

Angiographic & Cardiac System

CGO-3000



Product Introduction

1. HF 100kw Generator Control System(INDICO 100SP 100kW CPI)
2. Ergonomic Mechanical Design
3. 12" Image Intensifier(TH59432)
4. Excellent Image Quality
5. Three Acquisition Modes
6. Medical Monitor

Main Function

1. HF 100kw Generator Control System (INDICO 100SP 100kW CPI)
2. X-ray Tube Assembly (1.5 MHU)
3. Multi-Field Image Intensifier (TH59432)
4. Imaging System
5. Suspended Ceiling C-Arm
6. Cath-Table
7. Monitor 15"(1 set)
8. Monitor 19" (2 sets)

CGO-2100 FPD



Product Introduction

1. HF 100kw Generator Control System(INDICO 100SP 100kW CPI)
2. Ergonomic Mechanical Design
3. Paxscan 4030CB FPD
4. Excellent Image Quality
5. Three Acquisition Modes
6. Medical Monitor
7. Silent High-speed X-ray Tube with Oil Cooling System

Main Function

1. High Frequency Generator (100kW)
2. X-ray Tube Assembly (1500kHu)
3. Paxscan 4030CB FPD
4. Imaging System
5. 3-axes floor stand C-arm
6. Cath Table
7. Monitor 15"(1 set)
8. Monitor 19" (2 sets)

CGO-2100C



Product Introduction

1. HF 100kw Generator Control System(INDICO 100SP 100kW CPI)
2. Ergonomic Mechanical Design
3. 12" Image Intensifier(TH59432)
4. Excellent Image Quality
5. Three Acquisition Modes
6. Medical Monitor
7. Silent High-speed X-ray Tube with Oil Cooling System

Main Function

1. H.F. 100kw generator (touch screen control)
2. High-speed starter
3. X-ray Tube Assembly VARIAN G-1582(1500KHU)
4. Multi field image intensifier(TH59432)
5. DSA system
6. Cath table
7. C-arm
8. Monitor 15"(1 set)
9. Monitor 19"(2 set)
10. Monitor holder
11. Accessory