

Conventional Radiography X-ray Unit

F99-ICT



Product Introduction

1. Automatic and manual control programs: Programmable functions, provides 90 memory programs for radiography. The parameters can also be set manually.
2. X-ray tube protection: X-ray tube heat capacity and load capacity are managed and protected by computer programs, to allow maximum utilization of the tube capacity and radiographic parameters.
3. Anatomical exposure selection: Exposure parameters can be set based on patient body size and exposure positions. The parameters can also be corrected and stored manually.
4. Safety function: The error self-diagnosis function can determine problems quickly and accurately. It can execute protection automatically and display error code

Main Function

1. 500mA generator with microprocessor control
2. 4-way-floating table
3. APR with 90 memory programs