


Diagnostic Table

ZC15XY-1




Product Introduction

1. Table top Dimension: 2044×626mm, Moving range: Longitudinal 30cm extendable at the head end.
2. Tilting range: -15°~+90°
3. SID: 610mm ~ 910mm

Main Function

1. Table body
2. Control cabinet
3. Remote controller
4. Spot film device
5. Footrest

ZC25SY



Product Introduction

1. Control: Remote control
2. Table tilting range: +90°~-25°Speed : 4.4°/s
3. Length: 2100mm, Width: 625mm
4. Movement of spot film device: 740±10mm(0°~90°), 340±10mm(0°~25°)
5. Transfer movement of tabletop: 220±10mm
6. Footrest: Can be removed for gynecology examination
7. Collimator: Manual & Auto
8. Grid: r10, f120cm

Main Function

1. Table body
2. Control cabinet
3. Remote controller
4. Spot film device
5. Collimator
6. Footrest
7. Accessories

ZC35SY-1



Product Introduction

1. Dimension of table top: 2350×796mm
2. Vertical Movement of the Table: 822mm ~ 1262mm (when table is level)
3. Rotational Movement of the Table: +90°~ -35°
4. Focus and Film Distance (FFD): 1100mm, 1500mm (adjustable)
5. Table Movement (lateral): 250mm±10mm
6. SID: 1100MM, 1500MM
7. Table top-film distance: 77mm

Main Function

1. Table body
2. Control cabinet
3. Remote controller
4. Spot film device
5. Collimator
6. Footrest
7. Accessories