

High Frequency Mobile Unit

PX-100CLK



Main Function

1. High frequency 1.6kW generator
2. APR Function available
3. High image quality with low radiation exposure
4. Remote control available
5. CE Approved

XC30 H.F.



Main Function

1. High frequency 5kW generator
2. 3-fields, 9" Image Intensifier
3. Sensitive touch button, LCD display
4. CE approved

PX2000 H.F.



Main Function

1. High frequency 20kW generator
2. APR Function available
3. High image quality with low radiation exposure
4. Remote control available