

High frequency radiographic system

HF50-R50kw



Product Introduction

1. Normal radiography, bucky radiography, tilting radiography, Ion chamber AEC and APR mode can be selected.
2. Automatic heat capacity protection.
3. New-designed table, displays SID and rotation angle automatically.
4. High frequency design (25kHz).
5. Short exposure time (min 2ms).
6. High kV (150kVp), high mA (630mA), small size focus (0.6mm).
7. AEC+APR function.

Main Function

1. H.T. generator
2. H.T. cable 75kV/12m
3. X-ray tube assembly Toshiba E7843X
4. Radiographic table (including Ion Chamber)
5. X-ray tube stand
6. Collimator