

# Radiographic Table

## SC4-4



Table top	Dimension	2000×740×658mm
	Moving range	Longitudinal: 1200mm
		Transverse: 200mm
Lock mode	Electromagnetic lock	
Bucky	Moving range	Longitudinal: 500mm
	Lock mode	Electromagnetic lock
	Grid	Type: fixed/oscillating
		N40, fo100, r10
Film size	Max 356mm×432mm(14" X 17")	
Tube stand	Moving range	Longitudinal: 1200mm

## SC3-4



Table top	Dimension	2105×810mm
	Moving range	Longitudinal: 900mm
		Transverse: 240mm
Lock mode	Electromagnetic lock	
Elevating		488mm ~ 828mm
Bucky	Moving range	Longitudinal: 480mm
	Lock mode	Electromagnetic lock
	Grid	Type: fixed/oscillating
		N40, fo100, r10
Film size	Max 356mm×432mm(14" X 17")	
Ceiling suspended X-ray Tube	Moving range	Longitudinal: 2200mm
	Rotation	Vertical: 1200mm -45° ~ +90°
SID		1100mm
Tabletop-film distance		≤70mm

# Radiographic Table

## SC5-1



Table top	Dimension	2000×740×658mm
	Moving range	Longitudinal: 1200mm
		Transverse: 200mm
Lock mode	Electromagnetic lock	
Bucky	Moving range	Longitudinal: 500mm
	Lock mode	Electromagnetic lock
	Grid	Type: fixed/oscillating
		N40, fo100, r10
Film size	Max 356mm×432mm(14" X 17")	
Tube stand	Moving range	Longitudinal: 1200mm
	Rotation	Around lateral arm: +120° ~ -120°
		Around column axis: 0° ~ ±180°, can be fixed at ±90°, ±180°
SID		1100mm
Tabletop-film distance		≤60mm
Tomo program		8°, 20°, 40°
Tomo height		≥210mm